POWER POST

ENGLISH

BM4

BM4

PRESENTATION: 5 units

H 20mm C	DIMENSIONS:							
	Nº 0,5		Nº 1		Nº 2		Nº 3	
H	20,0mm	H	20,0mm	H	20,0mm	H	20,0mm	
B	1,4mm	B	1,6mm	B	1,8mm	B	2,0mm	
C	0,65mm	C	0,85mm	C	1,05mm	C	1,25mm	

DESCRIPTION: Power Post is a pin manufactured with a resinous die with a fiberglass base and epoxy resin which has a translucent presentation, being able to be used in regions which have greater esthetic requirements. The conic format of the pins allows a better adaptation to the radicular canals without requiring greater dental structure wear.

INDICATION: Power Post is indicated for reconstructing teeth which have undergone endodontic treatment and require restoration. It acts as intraradicular reinforcement of the dental structure and promotes retention for the final restorative material, or filling core, in the case of indirect restorations.

COMPOSITION: Fiberglass and epoxy resin.

PIN SIZE	
0,5/1	Molars and premolars; Upper and lower central and lateral incisors; Case of partial loss of dental structure; Reinforcement of nonvital teeth, weakened by excessive canal wear.
2/3	Cases with greater mechanical stress (significant loss of dental structure); Larger canals (with greater wear).

HOW TO USE: 1. Select the Pin in accordance with the canal diameter positioning the pin on the tooth radiograph; 2. Isolate the operating field; 3. Remove the gutta-percha from the canal using an appropriate instrument. Leave at least 3 to 5 mm of filling in the root apex, due to the presence of accessory canals, which must remain sealed to avoid bacterial contamination; 4. Prepare the root canal with suitable burs until adapting them to the format of the pin; 5. Take a periapical X-ray of the pin selected to check the preparation of the canal and condition of the filling. The pin must be perfectly adapted to the preparation executed, there must not be any spaces between the end of the pin and the gutta-percha; 6. If necessary cut the pin crown excess and put it in position again to check the adjustment; 7. Clean the pin with alcohol and dry it; 8. Apply the phosphoric acid gel 37% (Power Etching BM4) for 15 seconds inside the root canal and remaining dental structure. Wash thoroughly and remove the excess water with absorbent paper cones; 9. Apply a layer of Silane (Maquira) on the surface of all the pin and let it dry for a minute; 10. With a suitable bur, take the resinous cement to the canal; 11. Insert the pin in the canal and proceed with the photopolymerization of the cement (in accordance with the manufacturer's specification) by the surface and by means of the pin; 12. Note that the pin shall be inserted until the bottom of the prepared cavity and immobilized with accessory pins in the cases of larger canals.

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PREPARATION SEQUENCE WITH DRILL LARGO:

PRECAUTIONS: Wear safety goggles during the cutting of the pins and avoid cutting near the patient's face; Do not use scissors or pincers to cut the pin. It must be cut with diamond tipped burs in high rotation and under cooling; Use pincers to hold any fiberglass pin; Wear latex gloves to handle the product; The pins must be free of oily dirt, organic materials or any other contaminating agents at the moment of cementing to avoid failure of adhesion; Reprocessing forbidden.

WARNINGS AND/OR CONTRAINDICATIONS: It is not indicated in extensive oral rehabilitation and in patients with parafunctional habits; The technique of direct cores with fiber pins must only be applied in teeth with at least 2mm in height of coronary dentin.

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