PRESENTATION:
1 syringe with 2 g and 5 tips.

COLORS:
Yellow I Light I Clear I Yellow + I Light +.

INDICATION:
NATURES VENEER light-curing resin cement is indicated for the cementation of "dental contact lenses" and indirect cements that do not present opaque infrastructure. It provides permanent adhesion through the bonding and retention of the substrates through the photopolymerization system, being presented in a single component.

COMPOSITION:
Methacrylate monomers, HEMA, photoinitiators, preservative, stabilizer, fillers and pigments.

MODE OF USE:
Color proof: 1. Remove temporary veneer and temporary cement residue with brush/bowl and prophy paste. 2. Select a Natures Try-In color matched to the shade of the created veneer and adjacent teeth. 3. Dry the surface of the tooth to be restored, then place a small amount of Natures Try-In on the inside of the part and place it on the tooth. Any excess paste should be removed. 4. Evaluate the color matching effect created and if necessary try other colors, following the recommendation above. 5. Remove the piece carefully and wash it and the tooth surface thoroughly.

Preparation of the prosthetic piece: 1. The inner surface of the veneer can be blasted (e.g. aluminum oxide) to increase micro retention (laboratory procedure); for all-resin parts, this micro retention can be done with drills or aluminum oxide blasting. 2. For conditionable ceramics, apply hydrofluoric acid gel to its internal face, according to the time specified in the product's instruction manual. Rinse with plenty of water and dry. 3. With the internal surface clean, apply a layer of Silane for 10-15 seconds and allow to dry for 1 minute. 4. Apply the dental adhesive system as recommended in its instruction manual. Do not light cure.

Tooth preparation: 1. Ensure that the tooth surface and prosthetic piece are free of Natures Try-In and/or other substance. 2. Proceed to relative or absolute isolation according to professional preference and according to the teeth to be cemented. 3. Etch enamel and/or dentin surface for 15 seconds with 37% phosphoric acid gel. Wash and dry gently without dehydrating. 4. If the cervical endpoint is subgingival and relative isolation is preferred, perform a retraction maneuver with retractors. 5. Apply the dental adhesive system as recommended in its instruction manual. Do not light cure.

Cementation: After preparing the piece and tooth, proceed with the cementation itself. 1. When handling the cement, ensure that light from the reflector does not shine directly onto the cement to accelerate polymerization. 2. Place a thin layer of Natures Veneer onto the inner surface of the prepared veneer. 3. Position the piece on the prepared tooth and, while maintaining light pressure, remove excess cement with a micro-applicator or exploratory probe. 4. Prior to photoactivation, apply glycerin gel to the restoration margins to prevent an oxygen-inhibited layer. Light cure for 40 seconds at each margin (mesial, distal and cervical) and directly on the vestibular surface. **Note:** For photoactivation use halogen or LED equipment with a minimum power of 450 mW/cm² and blue light emission spectrum (400-500 nm). 5. Finishing and polishing can be performed with fine sandpaper discs, sandpaper strips, polishing pastes and felt discs.

WARNINGS AND PRECAUTIONS:
Exclusive use of the dental professional. Strictly follow the instructions for use. Avoid contact with eugenol-based products, as this may delay product hardening. Keep out of reach of children and pets. Do not use the product after its expiration date. Do not reuse the empty packaging. Irritation may occur with the dental pulp, if there is contact with it or in places close to it. In case of allergic reactions, discontinue use. In case of pulp exposure, protect with lining material. Avoid contact with the skin. In case of contact, remove the excess immediately with absorbent cotton soaked in water and rinse with plenty of water. Avoid contact with the eyes. If contact occurs, rinse with plenty of water and seek medical assistance immediately. It is recommended to adjust the reflector light so that it does not shine directly onto the restoration, promoting premature light curing of the cement. Use halogen or LED light equipment with a power equal to or greater than 450 mW/cm² and blue light emission spectrum (400 to 500 nm). Do not reuse the applicator tips to avoid cross infection. Verify the correct closure of the syringes after use. In case of allergic reactions to the product, discontinue use immediately.

STORAGE:
Keep the product away from direct light and heat. Temperature: 5 °C - 30 °C (41 °F - 86 °F).

DISPOSAL:
Leftovers or residues of used materials must not be reused, but must be disposed of in accordance with local legislation. After use, dispose of the tips according to local legislation.

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