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PRESENTATION:

Syringe with 4 g of the product, in the below listed colors:

Opaque: OA2 | OA3

Dentin: DA1 | DA2 dentin universal | DA3 | DA3,5 | DA4 | DB1 | DB2 | DB3 | DC2 | DC3 | D-Bleach

Enamel: EA1 | EA2 | EA3 | EA3,5 | EA4 | EB1 | EB2 | EB3 | EC2 | EC3 E-Bleach H | E-Bleach M | E-Bleach L

Translucent: T-Neutral | T-Blue | T-Yellow | T-Orange | T-Gray

DESCRIPTION:

Zirconfill is a nano-hybrid, light curable and radiopaque composite resin. Zirconfill has been developed to be used in the restoration of anterior and posterior teeth. It has a high content of inorganic filler (80% in weight). The filler content is mainly composed of zirconium/silica particles. The particle size distribution allows the adequate filling of the resin, contributing to its high mechanical resistance and wear resistance, and to the easy polishing; creating a restoration with smooth surface and high shine. It presents high workability and ideal consistence. It sets in the dental cavity, or in models, by the use of direct halogen light or LED.

TECHNICAL DATA:

Average particle size: 15.8 micrometers

Average nanoparticle size: 20 nanometers

Average percentage of filler in weight: 80%

INDICATIONS:

Zirconfill has been developed to be used in the restoration of anterior and posterior teeth. It is indicated to be used in: Direct restoration of the enamel and dentin, of all classes (I, II, III, IV e V); Esthetic corrections; Direct composite veneers; Tooth fragment bonding; Indirect restorations (inlays e onlays); Confection of core resin buildups; Anterior teeth stabilization.

COMPOSITION:

Mixture of Monomers (Bis-GMA, Bis-EMA, TEGDMA and UDMA), Photo initiator, Diatomite, Silica, Zirconia/Silica, Mixed Oxide and Colorants.

MODE OF USE:

1. Perform the prophylaxis of the teeth's surfaces. For that procedure, avoid the use of substances containing oil vehicle and/or fluoride. 2. Choose the appropriate color of the restorative material. For purposes of the color selection, the teeth must be clean and hydrated, without any isolation. 3. Perform the isolation of the operative field. The isolation must guarantee a moisture-free, clean and visible operative field, during all treatment. Preferably, use rubber dam isolation. 4. Prepare the cavity according to the recommendation for the execution of direct adhesive restorations. 5. Perform pulp protection, if necessary, with lining material based on calcium hydroxide and/or glass ionomer, according to the professional's clinical evaluation and, following the instructions for use of these materials. 6. For anatomical restoration of the teeth and, in order to minimize the excess resin to be removed, universal or pre-formed polyester matrixes must be used. The use of wedges will make the adaptation of matrixes and the achievement of proximal contact easier. 7. Perform the acid etching, following the instructions for use of the material that will be used. 8. Apply the adhesive system, according to the manufacturer's recommendations. 9. Perform the cavity restoration, using an adequate instrument to apply the increments of Zirconfill, carefully adapting them to the cavity walls. 10. We recommend the technique of insertion and polymerization by increment. This technique allows a more uniform polymerization. Use the maximum increment thickness and polymerization time, described in the chart below:

COLORS	Maximum increment thickness	Polymerization time
Dentin/Opaque	1.5 mm	30 s
Enamel/Translucent	2.0 mm	20 s

11. Light-cure each increment, as described in the chart above. Note that the presented values are ideal when light curing units with wavelength of 400 a 500nm are used. These values must be adjusted according to the equipment used by the professional. 12. If matrixes are used for the confection of the restoration, at the end of the execution and after removing the matrix, the last layer of resin must be submitted to an additional polymerization. 13. Check the occlusion with articulation paper. Carefully adjust the occlusion. 14. Perform finishing and contouring of the surface of the restoration. Finishing, contouring and oclusal adjust must be performed with instruments of choice of the professional, such as fine grit diamond burs, scalpel blades, silicone points, etc. Proximal surface contouring must be performed with sandpaper strips. 15. Perform polishing with sandpaper disc, felt discs and polishing pastes, of choice of the professional.

PRECAUTIONS:

The product is manufactured to be used only as exposed in the instructions for use. The use in a different way then described in the instructions for use is of single and exclusive responsibility of the professional. Product is for Professional dental use. These resins are light sensitive. For this reason, too intense ambient lights may cause the premature polymerization of the material. Avoid intense light in the operative field. After removing the resin from the syringe, use it fast. And immediately replace the syringe cap. After use, the product must be properly closed. Contact with saliva or blood during the application of resin might damage the restoration. The use of rubber dam or cotton rolls, along with saliva ejector, is recommended for adequate isolation. Avoid the use of lining or temporary materials based on eugenol, because they may interfere in the setting of the resin. Use masks, gloves and goggles for appropriate protection. It is also recommended that patients wear protective goggles too. Do not allow product to come contact with oral tissues or the skin. Do not use the product with expired validity. Cationic mouth solutions may cause discoloration, same as plaque disclosures and chlorhexidine. Do not use in patients with hypersensitivity to the components of the product. Keep out of reach of children. Single use product. Prohibited reprocessing.

RESTRICTIONS, WARNINGS OR CONTRAINDICATIONS:

There is no contraindication, when Zirconfill is used by a dental professional, according to the instructions for use. However, in patients with history of allergy to any components of the formula, the use is not recommended. In case of allergic reaction, stop the use and seek medical attention.

EXPIRATION DATE:

The expiration date is described in the packaging and represents the maximum deadline for use of the product. Do not use after the expiration date.

STORAGE AND TRANSPORT:

Product must be kept in the original packaging, which must be closed after use. Keep it protected from direct light incidence. Keep at temperatures between 5 and 30 °C, in a dry place. If product is kept under refrigeration, it must be removed from the refrigerator before use, to allow it to return to room temperature. Under adequate storage conditions, Zirconfill has a shelf life of 36 months.

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Zirconfill

MADE IN BRAZIL **maquira**
Dental Group

Fabricante | Manufacturer:
Maquira Indústria de Produtos Odontológicos S.A.
CNPJ: 05.823.205/0001-90
Endereço | Dirección | Address:
Av. Melvin Jones, 773
Pq. Ind. Bandeirantes
Maringá - PR | 87070-030 | Brasil
Resp. Técnico | Technical Resp.:
Cristiane Furlan | CRF-PR: 14.038
Reg. ANVISA N°: 80322400121

SAP:
Atendimento ao cliente
Atención al cliente
Customer service
0800 726 5848
sap@maquira.com.br